Advanced analytical methods optimize rAAV process development outcomes and improve yields and full particle ratios

H Martins¹, G Ganjam¹, F Leseigneur¹, S Pasupuleti¹, M Wachowius¹, E Ayuso¹

[1] Siegfried DINAMIQS AG, Department of Analytical Development, Schlieren, 8952, Switzerland

Contact: helena.martins@dinamiqs.com

Introduction

Recombinant adeno-associated viruses (rAAVs) are the leading vehicles for gene therapy, but large-scale manufacturing that balances high yield and quality remains a major challenge. A key difficulty is minimizing empty particles while maximizing full capsids containing the gene of interest. To address this, we employed mass photometry (MP) as a rapid and precise method to distinguish empty, partial, and full rAAVs based on their mass distribution^[1]. We optimized upstream production parameters and downstream purification using MP, identifying optimal transfection conditions that yielded the highest proportion of full particles and improving purification by selecting fractions enriched in full capsids. Nanopore sequencing^[2] further confirmed the identity and purity of rAAV preparations, revealing differences in genome heterogeneity. Our findings demonstrate the critical role of advanced analytics in optimizing rAAV production and enhancing gene therapy manufacturing.





Mass Photometry Samux^{MP} (Refeyn)

Nanopore Sequencing Mk1C (Oxford Nanopore Technologies)

DINAMIQS

a Siegfried company

Results and Discussion

MP evaluates **USP** parameters

To optimize the genome titer and full-to-empty capsid ratio in rAAV5 production, we investigated key upstream parameters in a triple transfection process using HEK293 cells grown in suspension. Cells were transfected with the same plasmids for the gene of interest (pGOI) and the rAAV5 capsid (pRep2Cap5), but either pHelper A or pHelper B was used. Two stoichiometries were tested: a 1:1:1 and a 2:1:0.5 ratio of pRep2Cap5:pGOI:pHelper. Following an optimized clarification process, vector genome titers were quantified by ddPCR, and the percentage of full rAAV particles was determined using mass photometry (SamuxMP).

Table 1: Overview of the experimental conditions, vector genome titers and % full particles in the clarified harvests.

Condition No.	Plasmid ratio (pRep2Cap5:pGOI:pHelper)	pHelper	% Full	Titer ddPCR (vg/ml)	Total VGs
1	1:1:1	pHelper A	8.8%	2.49e+11	7.47e+12
2	2:1:0.5	pHelper A	8.4%	3.39e+11	1.02e+13
3	1:1:1	pHelper B	16.5%	2.18e+11	6.54e+12
4	2:1:0.5	pHelper B	8.2%	2.02e+11	6.06e+12

1:1:1 Helper A			2:1:0.5 Helper A			
200 -		8.8% Full				8.4% Full

MP evaluates **DSP** parameters

We aimed to develop an efficient purification process for rAAV8, focusing on the separation of empty and full capsids using anion exchange chromatography (AEX). Mass photometry (MP) was employed to monitor the percentage of full capsids and detect potential impurities during the step-gradient elution process.

Empty and full rAAV8 preparations by CsCl gradient ultracentrifugation



Figure 2: Efficient removal of empty capsids from an rAAV8 preparation measured by MP.



Figure 3: AEX-process chromatogram. Load material (LM) and four elution fractions

Eluate 1

-0.03436 σ 0.00141 954 counts (90%) Skewness: 0.000

-0.03608 σ 0.00094 364 counts (50%) Skewness: 0.000

-0.03

Figure 4: Mass histograms of

highlighted in orange.

collected fractions with % full particles

0.00

-0.025

0.000

0.2% Full

Empty: 975 (99.8%) Full: 2 (0.2%) Ambiguous: 0 (0.0%)

-0.075

28.1% Full

-0.050

-0.04978 o 0.00098 counts (31%) ewness: 0.000

Empty: 406 (62.3% Full: 246 (37.7%)

-0.09

-0.06

Ratiometric contrast

Ratiometric contrast

Eluate 2

Nanopore Sequencing of Final Products

To investigate thoroughly genome heterogeneity, Nanopore sequencing was employed to analyse the identity of two rAAV preparations. Both contained the same GOI but were produced with different pHelper and pRepCap. NGS was also used to detect other residual nucleic acids.





Figure 1: Mass histograms of four clarified harvest samples from cells transfected with different pHelper and two different plasmid stochiometries. The % full particles is highlighted in orange.

The mass histograms (Figure 1) showed that cells transfected with pHelper B at a 1:1:1 ratio produced the highest percentage of full rAAV particles (16.5%) compared to 8.2% with the 2:1:0.5 ratio. Although pHelper A resulted in higher overall productivity (Table 1), the percentage of full capsids averaged 8% for both stoichiometries.

(E1-E4) were collected and characterized by MP.

Empty: 1042 (81.5%) Full: 237 (18.5%) Ambiguous: 0 (0.0%)

-0.075

Empty: 294 (71.9% Full: 115 (28.1%)

-0.075

Empty: 346 (88.9%) Full: 43 (11.1%) Ambiguous: 0 (0.0%

-0.075

31.0% Full

Editor a a a a a a

11.1% Full

-0.100

-0.100

Load material

-0.03584 σ 0.00132 994 counts (71%) Skewness: 0.000

-0.025

-0.03631 σ 0.00080 278 counts (5 Skewness: 0.0

-0.025

-0.03621 σ 0.00106 317 counts (43%) Skewness: 0.000

-0.025

-0.050

-0.050

-0.050

Ratiometric contrast

Ratiometric contrast

Eluate 4

Ratiometric contrast

Eluate 3

0.000

0.000

18.5% Full



Figure 5: A) Simplified sequencing workflow. B) Coverage analysis of reads aligning to the vector genome from ITR-to-ITR. C) Coverage analysis of reads aligning to the full pGOI. The percentage of mapped reads aligning to the backbone (red) and the vector genome (gray) is displayed. D) Ratios of reads mapped to the respective production plasmids or unmapped reads.

Our in-house bioinformatics pipeline revealed that Vector A showed significantly lower DNA contamination than Vector-B. Specifically, rAAV-A had 1.8% from backbone packaging, 1.2% from pRepCap, 0.4% from pHelper, and 5.5% from host cell

the final elution (E4), however, a notable proportion of lowermass impurities was also present.

Among the elution fractions, E2 and E3 exhibited the highest

enrichment of full particles, with 31.0% and 28.1%, respectively.

A small population of full capsids (11.1%) was also detected in

DNA, whereas rAAV-B had 8.4%, 4.8%, 1.7%, and 19.7%, respectively.

Conclusions

Our findings underscore the importance of identifying the right analytical tools to understand and optimize each step of rAAV manufacturing and highlight how an agile integration of advanced analytics accelerates process development:

- We effectively used MP to distinguish between empty, partial, and full rAAV particles to identify optimal upstream production parameters.
- MP enabled rapid and accurate monitoring of full-to-empty particle ratios during downstream purification via AEX. This facilitated the selection of the best elution fractions while identifying impurities in other fractions.
- Nanopore sequencing provided detailed insight into genome identity and residual levels from other nucleic acids.

References

1. Wagner C., et al. Quantification of Empty, Partially Filled and Full Adeno-Associated Virus Vectors Using Mass Photometry. Int. J. Mol. Sci. 2023, 24, 11033, doi: 10.3390/ijms241311033. 2. Marco T Radukic, et al. Nanopore sequencing of native adeno-associated virus (AAV) singlestranded DNA using a transposase-based rapid protocol, NAR Genomics and Bioinformatics, Volume 2, Issue 4, December 2020, Iqaa074, https://doi.org/10.1093/nargab/Iqaa074

Acknowledgements

1. This work is funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or HaDEA. Neither the European Union nor the granting authority can be held responsible for them. 2. This work is funded by the UK Research and Innovation (UKRI) under the UK government's Horizon Europe funding guarantee grant numbers 10080927 (The Francis Crick Institute Limited), 10079726 (King's College London), 10082354 (Muscular Dystrophy Group of Great Britain and Northern Ireland) and 10078461 (University College London). 3. This work has received funding from the Swiss state Secretariat for Education, Research and Innovation (SERI).

